

ex
extriCARE®



extriCARE® 2400 (MI397)
Negative Pressure Wound Therapy System
Distributor Sell Sheet

extriCARE® 2400

Negative Pressure Wound Therapy System



The extriCARE 2400® pump is a lightweight, portable, battery-powered negative pressure wound therapy pump.

This device, used in conjunction with our anatomically fitted line of wound dressings, may promote wound exudates, infectious material, and tissue debris from the wound bed using continuous and/or intermittent negative pressure.

Features and Benefits:

- Effective:** Powerful negative pressure, may promote faster wound healing
- Reliable:** Rechargeable battery, 30 hours of continuous battery life with a full charge
- Safe:** Airtight, disposable canisters that come in two sizes, gelling agent solidifies wound exudates

Device Specifications:

Product Code:	<i>EC2400 (MI397)</i>
Weight:	8.6 oz.
Battery Type:	Lithium (rechargeable)
Vacuum Modes:	Continuous or intermittent
Charging Time:	< 3.5 hours
Dimensions:	Length: 3.35" (8.5 cm) Height: 5.67" (14.4 cm) Width: 1.46" (3.7 cm)
Input:	100-240 V 50/60 Hz 0.4 A